

Mid-fire body clays.

	Clay	Filter	Bisque (Orton cones)	Glaze (Orton cones)	Characteristics
MFQ	Mid-fired quartz	120#	06 to 05	1 to 3	A fine, white mid-firing clay. It is an excellent clay for throwing, pressing, jigger jollyng, and so on. MFQ offers stoneware performance when fired in the range of Orton cones 1 to 5. Excellent glaze fit and underglaze intensity.
SWE		120#	6	02 to 6	A versatile cream to buff, earthenware/stoneware clay. It is very tolerant to firing conditions. Glaze fit at 1200°C. It offers fine texture and good plasticity. Great for fine detail and excellent plasticity.
LBCI	Cool Ice	120#	07 to 05	4 to 6	A truly magnificent, translucent, low firing, white porcelain. The firing range is 1186 to 1222. For best whiteness and translucency fire in reduction.

Notes:-

1. Unless otherwise specified the bodies are filter pressed, vacuum de-aired and packed as bags of approximately 10 kg weight, i.e. 100 bags to the tonne.
2. Conditions and production techniques vary widely from pottery to pottery; consequently, we strongly advise thorough evaluation and testing of any clay body to ensure its suitability.
3. Clayworks reserves the right to vary body specifications and composition, without notification.
4. The information provided is the latest available and supersedes all earlier publications.